


Fast and Easy Detection of FII

Risk Factors of Genetic Thrombosis:

An *in vitro* diagnostic test that allows the detection and genotyping of a single point mutation (G20210A) of the human Factor II gene, referred to as Factor II mutation, directly from human whole blood.




Reliable genotyping of FII mutation:

LC-FII-LP-24™ kit improves dramatically the efficiency of the process.

- ➔ **RAPIDITY:** 1 minute hand-time & 30 minutes amplification + detection
- ➔ **SIMPLICITY:** No laborious DNA extraction procedure
- ➔ **EFFICIENCY:** Perfect distinction of the heterozygote

WORKFLOW



PLATFORM

➔ LC-Genie III™



KIT COMPOSITION

LC-FII-LP-24™	Positive Control Negative Control LaCAR reaction Mix 1 LaCAR reaction Mix 2	24 PCR reactions
----------------------	--	------------------